




Interventional Radiology

Free Shipping on Orders \$99+ *



Phone (877) 220-4567 - Fax (786) 231-0336
www.edm-us.com - info@edm-us.com

General Purpose Ultrasound Probe Covers

Latex-free Probe Cover Kit by EDM

- * Best value for ultrasound-guided interventions
- * Sterile, single-use
- * Latex-free polyethylene
- * Telescopic fold

Item	Description	Pckg	w/ Sterile Gel	Price
101.0016	5.5" x 12"	25/box	Yes	\$99
101.0017	5.5" x 32"	25/box	Yes	\$109
101.0015	5.5" x 48"	25/box	Yes	\$119
101.0019	5.5" x 12"	100/box	No	\$270
101.0018	7" x 48"	50/box	No	\$208



Best Value

Recommended for Vascular Access (IV or catheter placement), pain management (injections, PRP, Stem cells), and operating room (nerve blocks, regional anesthesia)

Probe Cover Kit by Protek

- * Sterile, single-use
- * Telescopic fold
- * Latex-free



The Soft material is made of latex-free polyurethane. The Clear material is made of latex-free polyethylene.

PullUp Probe Covers

- * PullUp saves time when installing the cover
- * Kit with sterile gel & elastic bands
- * Sterile, single-use
- * Latex-free



Item	Description	Pckg	w/ Sterile Gel	Price
101.0005	5.5" x 36" - Soft material	24/box	Yes	\$161
101.0006	5.5" x 48" - Soft material	24/box	Yes	\$139
101.0007	6" x 68" - Soft material	24/box	Yes	\$243
101.0013	6" x 48" - Clear material	100/box	No	\$249
101.0014	6" x 48" - Clear material	24/box	Yes	\$139

Item	Description	Pckg	Price
101.0001	5" x 24" - Soft material - #2099	24/box	\$169
101.0002	5" x 36" - Soft material - #2095	24/box	\$191
101.0003	5" x 48" - Soft material - #2097	24/box	\$179
101.0004	5" x 48" - Clear material - #1096	24/box	\$142

Latex-free Probe Covers by Microtek

- * Stretchable latex-free material. Fits small and large transducers
- * Sterile, single-use
- * Kit includes latex cover, sterile gel and elastic bands



Request a free Sample Today

Item	Description	Pckg	w/ Sterile Gel	Price
101.0031	5.5 x 48"	20/box	Yes	\$139
101.0032	5.5 x 24"	20/box	Yes	\$159

SaferSonic Probe Covers

- * Covers with pre-gelled window. Does not require gel inside!
- * Easy clean up: no gel is left on the transducer
- * Sterile, single-use
- * Telescopic fold
- * Latex-free



Item	Description	Pckg	w/ Sterile Gel	Price
101.0020	7" x 12"	50/box	No	\$279
101.0021	7" x 12"	50/box	Yes	\$349
101.0022	7" x 47"	50/box	No	\$319
101.0023	7" x 47"	50/box	Yes	\$369
101.0024	7" x 98"	50/box	No	\$349

Surgical Ultrasound Probe Covers

Accordion-folded Surgical Covers

- * Extended length suitable for OR
- * Sterile, single-use
- * Latex-free
- * Cover made of polyethylene
- * Boot made of polyurethane



Item	Description	Pckg	Price
103.0001	6" tapered to 3" x 96" w/ boot, gel and bands	12/box	\$151
103.0002	6" tapered to 3" x 96" w/ boot and bands	12/box	\$139
103.0003	6" tapered to 3" x 96" w/ bands	12/box	\$129

Two-Piece Surgical Covers

- * Extended length suitable for OR
- * Sterile, single-use
- * Latex-free
- * Cover made of polyethylene
- * Boot made of polyurethane
- * 2 pieces must be assembled before procedure



Item	Description	Pckg	Price
103.0009	5" tapered to 3" x 96" Pull-Up w/ Soft Boot	24/box	\$239
103.0010	6" tapered to 3" x 96" Telescopic w/ Soft Boot	24/box	\$239

Pull-Up Surgical Covers

- * Two parts factory-sealed together:
Patient end is soft polyurethane (48")
Equipment end is clear polyethylene (48")
- * Extended length suitable for OR
- * Sterile, single-use
- * Latex-free



Item	Description	Pckg	Price
103.0008	Pull-Up Cover 5" x 96" w/ gel	24/box	\$255

Telescopic Surgical Covers

- * Extended length suitable for OR
- * Sterile, single-use
- * Latex-free
- * Cover made of clear polyethylene



Item	Description	Pckg	Price
103.0007	Telescopic Covers 5.5" x 96" w/ gel	24/box	\$215

Sterile Tablet Covers

- * Protects your tablet in the OR
- * Peelable adhesive tape on the back
- * Sterile, single-use
- * Latex-free



Item	Description	Pckg	Price
107.0001	Tablet Covers 8" x 10"	24/box	\$70

TEE Probe Covers

- * For transesophageal echocardiography
- * Protects patient-monitoring equipment during surgery
- * Cover rolled over plastic insert
- * Single-use
- * Latex-free



Item	Description	Pckg	Sterile	Price
104.0001	Kit 1.75" x 48" w/ bite block, gel, and tape	12/box	Yes	\$213
104.0002	Kit 1.75" x 48" w/ Tape	12/box	No	\$168
104.0003	Kit 1.75" x 48" w/ bite block, gel, and tape	12/box	No	\$188

Portable Ultrasound System Covers

- * Designed for portable ultrasounds
- * Peelable adhesive tape on the back
- * Sterile, single-use
- * Latex-free



Item	Description	Pckg	Price
107.0002	Portable System Covers 28" x 17.5"	24/box	\$149

Biopsy Gun Covers

- * Designed for reusable biopsy guns
- * Protects equipment, patients and staff from contaminants
- * Pull-Up system for easy insertion
- * Sterile, single-use
- * Latex-free




Item	Description	Pckg	Price
107.0003	Biopsy Gun Covers 3.75" x 17"	24/box	\$175

Endocavity Ultrasound Probe Covers

Latex Probe Covers - Non-Sterile

- * Ideal for routine transvaginal or transrectal examinations
- * Rolled
- * Individually wrapped (Except for #102.0009)

Item	Description	Pckg	Price
102.0006	0.8" x 7.9"	24/box	\$64
102.0007	1" x 7.9"	24/box	\$64
102.0008	1.5" x 7.9"	40/box	\$45
102.0009	1.5" x 7.9" - BULK	100/box	\$35



Latex Probe Covers - Sterile



- * Ideal for routine transvaginal or transrectal interventions
- * Rolled
- * Individually wrapped

Item	Description	Pckg	Price
102.0004	0.8" x 7.9"	24/box	\$70
102.0018	1.5" x 7.9"	20/box	\$35
102.0019	1.5" x 7.9"	50/box	\$75
102.0020	1.5" x 7.9"	100/box	\$125



Parker Eclipse Latex-free Probe Covers

- * Pre-Gelled Probe Covers
- * Non-sterile
- * Latex-free polyethylene

Item	Description	Pckg	Price
102.0014	Eclipse 2.37" tapered to 1.75" x 9.5"	100/box	\$62
102.0017	Eclipse 3D 2.52" tapered to 1.75" x 9.5"	100/box	\$85



Microtek PIP Probe Covers

- * Ideal for latex-sensitive patients and staff-members
- * Latex-free Polyisoprene
- * Rolled
- * Individually wrapped

Item	Description	Pckg	Sterile	Price
102.0022	1.4" x 12" (3.6 x 30 cm)	10/box	Yes	\$45
102.0023	1" x 12" (2.6 x 30 cm)	100/box	No	\$155
102.0024	1" x 12" (2.6 x 30 cm)	10/box	Yes	\$45
102.0025	1.2" x 12" (3 x 30 cm)	10/box	Yes	\$48



SlimFit Probe Covers with PullUp

- * Suitable for prolonged exams
- * Latex-free polyurethane
- * Telescopic with PullUp

Item	Description	Pckg	Sterile	Price
102.0021	4" tapered to 1.25" x 12" - No Gel	24/box	No	\$74
102.0010	4" tapered to 1.25" x 12" - Kit with Sterile Gel	24/box	Yes	\$188
102.0011	4" tapered to 1.25" x 24" - Kit with Sterile Gel	24/box	Yes	\$225
102.0026	4" tapered to 1.25" x 12" - No Gel	200/box	No	\$580



Latex-free Probe Covers

- * Ideal for transvaginal or transrectal exams of allergic patients
- * Flat fold
- * Latex-free polyurethane

Item	Description	Pckg	Sterile	Price
102.0001	5" tapered to 2" x 24" - Kit with Sterile Gel	24/box	Yes	\$158
102.0016	4" tapered to 1.25" x 12" - No Gel	10/box	No	\$26
102.0003	4" tapered to 1.25" x 12" - No Gel	50/box	No	\$88
102.0002	5" tapered to 2" x 24" - No Gel	24/box	No	\$59



Do you store your transducer in accordance with The Joint Commission guidelines?

"Store the device in a manner that will protect from damage or contamination and that is consistent with national guidelines and manufacturers' recommendations [...]"

Probe Storage Covers

- * Keeps your transducer from contaminants when not in use
- * Non-Sterile
- * Latex-free
- * Fits all endocavity transducers
- * Features a designated space for applying disinfection label
- * Not intended for ultrasound scanning

Item	Description	Pckg	Price
105.0001	4" x 16" - Flat fold	400/box	\$369
105.0002	4" x 16" with PullUp	200/box	\$369




**MORE
INFORMATION
HERE**



Latex Probe Cover Kit with Gel

- * Ideal for procedures where sterility is required
- * Rolled
- * Kit contains an individually wrapped cover and a gel packet
- * Extremely stretchable

Item	Description	Pckg	Sterile	Price
101.0028	1.5 x 7.9" with Gel	40/box	Yes	\$147



Ultrasound Gels

Krystal Ultrasound Gel by EDM

- * Non-sterile, bacteriostatic
- * High-viscosity
- * Dye-free, paraben-free, silicon-free



Item	Description	Pckg	Price
401.0015.12	8.5 oz (250 ml) dispenser bottle	12/box	\$18
401.0015.24	8.5 oz (250 ml) dispenser bottle	24/box	\$36
401.0015.40	8.5 oz (250 ml) dispenser bottle	40/box	\$60

Parker Aquasonic 100

- * Non-sterile
- * High-viscosity
- * Blue-colored

Acoustically correct
for the broad range of
frequencies used



Item	Description	Pckg	Price
401.0001	60g Tube	12/box	\$20
401.0002	8.5 oz (250 ml) dispenser bottle	12/box	\$24
401.0003	1 l dispenser bottle	6/box	\$54
401.0004	5 l SonicPack	4/box	\$84
401.0013	20 ml packet - Non-sterile	100/box	\$39

Parker Scan

- * Non-sterile
- * Low-viscosity
- * Green-colored

For medical ultrasound
procedures when a less
viscous gel is desired



Item	Description	Pckg	Price
401.0010	8.5 oz (250 ml) dispenser bottle	12/box	\$45
401.0011	1 gallon bottle	ea	\$28
401.0012	ScanPac w/ 4 gallons, 1 pump and 2 bottles	ea	\$99

Thermasonic Gel Warmer

- * Warms gel to body temperature
- * "Thermal cutoff" provides protection from overheating
- * Can be wall-mounted or used on a countertop

Item	Description	Pckg	Price
403.0001	Single-bottle	ea	\$125
403.0002	3-bottle	ea	\$200
403.0003	3-bottle w/ LCD temperature	ea	\$225



EcoVue Ultrasound Gel

- * Non-sterile, bacteriostatic
- * High-viscosity
- * Dye-free, paraben-free, silicon-free



Item	Description	Pckg	Price
401.0019	8.5 oz dispenser bottle	12/box	\$28

Parker Aquasonic Clear

- * Non-sterile
- * High-viscosity
- * Dye-free

Acoustically correct
for the broad range of
frequencies used



Item	Description	Pckg	Price
401.0006	8.5 oz (250 ml) dispenser bottle	12/box	\$24
401.0007	1 l dispenser bottle	6/box	\$54
401.0008	5 l SonicPack	4/box	\$84
401.0009	60 g Tube	12/box	\$20

Sterile Ultrasound Gel

- * Sterile, single-use
- * High-viscosity

Recommended for
interventional or
endocavity procedures



Item	Description	Pckg	Price
402.0002	EDM Sterile gel - 20 ml	48/box	\$55
402.0001	Parker Aquasonic 100 - 20 ml	48/box	\$65

Femina Lubricating Gel by EDM

- * Non-sterile
- * Hypoallergenic, non-greasy
- * Water-based & water-soluble
- * Alcohol-free, salt-free, scentless



Item	Description	Pckg	Price
407.0001	250 ml	25/box	\$75

Lubricating Gel by EDM



- * Hypoallergenic, non-greasy
- * Water-based & water-soluble
- * Alcohol-free, salt-free, scentless



Item	Description	Sterile	Pckg	Price
407.0002	5 g packet	No	100/box	\$15
407.0003	5 g packet	Yes	100/box	\$20
407.0004	42 g tube	No	6/box	\$9
407.0005	82 g tube	No	12/box	\$19

HR Sterile Lubricating Jelly

- * Sterile
- * Water-soluble
- * Hypoallergenic
- * Fragrance-free, salt-free, alcohol-free



Item	Description	Pckg	Price
405.0001	4 oz tube	12/box	\$21
405.0002	2 oz tube	12/box	\$19.2
405.0003	1.25 oz packet	48/box	\$48.8
405.0007	5 g packet	144/box	\$13.8

Surgilube Sterile Lubricant

- * Sterile
- * Water-soluble
- * Hypoallergenic, antimicrobial, bacteriostatic
- * Fragrance-free, salt-free, alcohol-free



Item	Description	Pckg	Price
406.0002	2 oz tube	12/box	\$28.5
406.0004	4.25 oz tube	12/box	\$43
406.0006	5 g foilpac	144/box	\$34.9

Parker Protex Wipes & Spray

- * Disinfects in 4 minutes
- * Alcohol, bleach, solvent-free
- * Can be used on surfaces
- * Approved by major OEM



Item	Description	Pckg	Price
601.0009	Canister - 100 wipes	6/box	\$79
601.0007	Softpack - 60 wipes	12/box	\$65
601.0004	12 oz bottle	12/box	\$84
601.0005	32 oz trigger-bottle	6/box	\$65

Sono Ultrasound Disinfection Wipes

- * Kills 99.9% of bacteria and viruses in 15 seconds
- * Alcohol, bleach, solvent-free
- * Wide spectrum: kills MRSA, HIV, HPV and staph
- * Approved by major OEM




Item	Description	Pckg	Price
601.0001	Softpack - 50 wipes	12/box	\$75
601.0003	Canister - 80 wipes	6/box	\$65

Disinfection Soaking Station

- * For high-level disinfection of ultrasound transducers
- * Can be wall-mounted or placed on a countertop
- * Dimensions: 11" x 7" x 7" (excluding the stand)

Item	Description	Pckg	Price
602.0001	Soak cup w/ stand	ea	\$89
602.0002	Replacement cups	3/box	\$45



High-level disinfection of ultrasound transducers

Accepts Glutaraldehyde, OPA and hydrogen peroxide

Surgical Equipment Drapes

Surgical Camera Drapes


- * Designed for endoscopy, arthroscopy
- * Universal design can fit all the cameras
- * Adhesive tape to secure the cover on the endoscope
- * Latex-free polyethylene
- * Sterile, single-use



Item	Description	Pckg	Price
108.0001	Camera drape W/ insertion ring - 6" x 98"	20/box	\$60
108.0002	Camera drape W/ telescopic fold - 6" x 98"	20/box	\$60
108.0003	Camera drape W/ accordion fold - 6" x 98"	20/box	\$60

Rectangular Banded Bags


- * Ideal to protect all types of surgical equipments
- * Fits a wide variety of OR equipments
- * Sterile, single-use
- * Latex-free polyethylene



Item	Description	Pckg	Price
110.0001	Rectangular Banded Bag 32" x 40"	100/box	\$215
110.0002	Rectangular Banded Bag 35" x 60"	85/box	\$215
110.0002	Rectangular Banded Bag 40" x 48"	90/box	\$225

Circular Dome Bags

- * Ideal to protect all types of surgical equipments
- * Fits a wide variety of OR equipments
- * Sterile, single-use
- * Latex-free polyethylene



Item	Description	Pckg	Price
109.0001	Circular Dome Bag 35" x 18"	100/box	\$195
109.0002	Circular Dome Bag 40" x 20"	100/box	\$195
109.0003	Circular Dome Bag 48" x 24"	100/box	\$215
109.0004	Circular Dome Bag 52" x 26"	100/box	\$215

Mini C-Arm Covers for Fluoroscan

- * Designed for Fluoroscan: covers the tube, the intensifier and the arm
- * Protects equipment, patients and staff from contaminants
- * Sterile, single-use
- * Latex-free

Item	Description	Pckg	Price
112.0001	Mini C-arm cover - 54" x 78"	24/box	\$209

Universal Light Handle Covers

- * Universal cover that fits most light handles
- * Sterile, single-use
- * Latex-free



Item	Description	Pckg	Price
111.0001	Universal light handle cover	120/box	\$120



Sony Thermal Papers

Item	Description	Pckg	Price
501.0001	UPP-110S - A6 B&W Standard	10 rolls/box	\$75
501.0002	UPP-110HA - A6 B&W High-contrast	10 rolls/box	\$139
501.0003	UPP-110HD - A6 B&W High-density	10 rolls/box	\$139
501.0004	UPP-110HG - A6 B&W High-gloss	10 rolls/box	\$139
501.0005	UPP-210SE - A4 B&W Standard	5 rolls/box	\$145
501.0006	UPP-210HD - A4 B&W High-density	5 rolls/box	\$180
501.0011	UPC-21S - Color paper - small	240 sheets/box	\$109
501.0012	UPC-21L - Color paper - Large	200 sheets/box	\$119
501.0013	UPP-84HG - High-gloss paper	10 rolls/box	\$109
501.0014	UPT-210BL - Blue transparent film	5 rolls/box	\$679

Mitsubishi Thermal Papers

Item	Description	Pckg	Price
501.0009	K65HM - High-density paper	4 rolls/box	\$49
501.0008	K61B - Standard paper	4 rolls/box	\$27
501.0010	K91HG - High-gloss paper	4 rolls/box	\$54



Sony UPX-898MD Hybrid Printer

- * Digital & analog printer with USB connection
- * Black & White A6 paper
- * Accepts UPP-110S, UPP-110HA, UPP-110HD & UPP-110HG



Item	Description	Pckg	Price
502.0001	UPX-898MD	ea	\$600

Sony UP-25MD Color Printer


- * Ideal for Medical Carts
- * Color A6 paper
- * Accepts UPC-21S & UPC-21L



Item	Description	Pckg	Price
502.0002	UP-25MD	ea	\$2,400

Sony Surgical Monitor Display

- * High resolution Full HD or 4K images
- * Easy, intuitive operation
- * Wide range of digital & analog signal inputs
- * Resists moisture and water
- * Easy clean edge-to-edge flush design



Item	Description	Pckg	Price
502.0003	Sony LMD-2435MD 24" Full HD	ea	\$2,555
502.0004	Sony LMD-2735MD 27" Full HD	ea	\$2,750
502.0005	Sony LMD-X2705MD 27" 4K	ea	\$4,950

Sony Medical Recorder

- * High Definition video recording
- * Simultaneous recording on internal HDD and external media
- * Burn files to DVD
- * Wide range of interfaces and network data transmission



Item	Description	Pckg	Price
502.0006	HVO-550MD Medical Recorder HD	ea	\$3,550
502.0007	HVO-550MD/FHD Medical Recorder Full HD	ea	\$4,125

Epson PP100 Disc Producer

- * Produces CDs, DVDs, and Blu-Ray discs
- * Holds 100 discs
- * Prints up to 60 discs an hour
- * Burns up to 30 CDs, 15 DVDs, or 8 Blu-rays per hour
- * 3 publishing modes



Item	Description	Pckg	Price
503.0001	Epson PP100N - Network capable	ea	\$4,850
503.0002	Epson PP100III	ea	\$4,850

Epson PP100 Ink & CDs



Item	Description	Pckg	Price
503.0003	Yellow ink cartridge for PP100	ea	\$42.95
503.0004	Black ink cartridge for PP100	ea	\$42.95
503.0005	Cyan ink cartridge for PP100	ea	\$42.95
503.0006	Light Cyan ink cartridge for PP100	ea	\$42.95
503.0007	Light Magenta ink cartridge for PP100	ea	\$42.95
503.0008	Magenta ink cartridge for PP100	ea	\$42.95
503.0009	White Inkjet Printable CD-R - 700MB	600/box	\$138

Biopsy Needle Guides & Training Phantoms

Ultrasound Needle Guides: Quicker and More Accurate Ultrasound-guided Procedures

General Purpose Needle Guides are tools designed to provide more precision during ultrasound-guided interventions. They improve needle visualization and patient safety. They also reduce procedure time and patient discomfort.

- * Reduces the number of attempts to reach a target
- * Procedures are less traumatic for patients
- * Improves needle placement accuracy
- * Provides a predictable needle path, with the help of on-screen guidelines

- * Sterile, single-use
- * Latex-free

Description	Pckg	Price
Director Guide Only	24/box	\$179
Director Procedure Kit	24/box	\$369
Retrofit Guide Only	24/box	\$159
Retrofit Procedure Kit	24/box	\$319



Radiology:

- * Biopsies
- * FNAs
- * Drainages

Pain Management:

- * Joint injections
- * Stem cell therapies
- * PRP injections

Anesthesia:

- * Nerve blocks
- * Injections
- * Infiltrations

Vascular Access:

- * PICC lines insertions
- * CVC placements

Endocavity Biopsy Needle Guides

- * Keeps the needle straight during prostate biopsies
- * Accepts 16 Ga & 18 Ga needles
- * Sterile, single-use
- * Latex-free



Item	Description	Pckg	Price
201.XXX1	Needle Guide Only	24/box	\$299
201.XXX2	Needle Guide Kit w/ 2 x 20 cm latex cover	24/box	\$359
201.XXX3	Needle Guide Kit w/ 2.6 x 20 cm latex cover	24/box	\$359
201.XXX4	Needle Guide Kit w/ 3.5 x 20 cm latex cover	24/box	\$359
201.XXX5	Needle Guide Kit w/ 4"-1.25" x 12" latex-free cover	24/box	\$379
201.XXX6	Needle Guide Kit w/ 5"-2" x 24" latex-free cover	24/box	\$379
201.XXX7	Needle Guide Kit w/ 2 x 20 cm latex cover & 5"-2" x 24" latex-free cover	24/box	\$416
201.XXX8	Needle Guide Kit w/ SlimFit PullUp latex-free cover & 4"-1.25" x 12" latex-free cover	24/box	\$439
201.XXX9	Needle Guide Kit w/ SlimFit PullUp latex-free cover	24/box	\$416

Needle guides are compatible with the following transducer brands: GE, Shimadzu, Samsung-Medison, Alpinion, Prosonic, Philips, Ultrasonix, Esaote, Siemens, Pie Medical, Terason, Sonoscape, BK Medical, Carestream, Hitachi-Aloka, Toshiba, Zonare, Mindray, International Biomedical, Chison, Supersonic Imagine. Please specify your transducer model upon order. Kit options vary upon probe model. Ask our customer service.

Training Phantoms

- * Safely improve your ultrasound and interventional skills
- * Material emulates human tissue properties under ultrasound: nerves, veins, arteries or lumps
- * Suitable for teaching, learning, training or assessing purposes
- * High-quality and enhanced durability

Item	Description	Pckg	Price
304.0001	Breast phantom w/ 11 lumps	ea	\$390
304.0002	Anesthesia phantom w/ 3 vessels & 3 nerves	ea	\$390
304.0003	Universal NAV phantom w/ 2 nerves, 2 arteries & 2 vessels	ea	\$390
304.0004	Thyroid phantom for thyroid biopsy training	ea	\$390
304.0006	Pediatric arm phantom	ea	\$420



MultiCore Automatic Biopsy Guns

- * Two firing modes, two triggers
- * Easy one-hand cocking
- * Sterile, single-use
- * Latex-free



Description	Pckg	Price
MultiCore Automatic Gun	15/box	\$329
SpringCut Semiautomatic Gun	15/box	\$275

MultiCore Automatic	Size	SpringCut Semiautomatic
211.0001	14 Ga x 11 cm	211.0017
211.0002	14 Ga x 15 cm	211.0018
211.0003	14 Ga x 20 cm	211.0019
211.0004	14 Ga x 25 cm	211.0020
211.0005	16 Ga x 11 cm	211.0021
211.0006	16 Ga x 15 cm	211.0022
211.0007	16 Ga x 20 cm	211.0023
211.0008	16 Ga x 25 cm	211.0024
211.0009	18 Ga x 11 cm	211.0025
211.0010	18 Ga x 15 cm	211.0026
211.0011	18 Ga x 20 cm	211.0027
211.0012	18 Ga x 25 cm	211.0028
211.0013	20 Ga x 11 cm	211.0029
211.0014	20 Ga x 15 cm	211.0030
211.0015	20 Ga x 20 cm	211.0031
211.0016	20 Ga x 25 cm	211.0032

SpringCut Semiautomatic Biopsy Guns

- * Stylet advance button remains locked during penetration
- * Pre-determined bioptic throw: 10-20 mm
- * Sterile, single-use
- * Latex-free



Chiba Needles

- * Easy-to-read centimetre marks and sliding stopper facilitate accurate depth placement
- * Quincke-type point facilitates and makes penetration less traumatic
- * Ultra-thin walled cannula
- * Sterile, single-use
- * Latex-free
- * 25/box
- * \$175



	18 Ga	20 Ga	21 Ga	22 Ga	23 Ga	25 Ga
5"	NA	NA	NA	NA	NA	212.0032
7"	212.0001	212.0008	212.0015	212.0021	212.0027	212.0033
9"	212.0002	212.0009	212.0016	212.0022	212.0028	212.0034
11"	212.0003	212.0010	212.0017	212.0023	212.0029	212.0035
15"	212.0004	212.0011	212.0018	212.0024	212.0030	NA
20"	212.0005	212.0012	212.0019	212.0025	212.0031	NA
25"	212.0006	212.0013	NA	NA	NA	NA
35"	212.0007	212.0014	212.0020	212.0026	NA	NA

Breast Localization Wire Needles

- * Easy-to-read centimetre marks facilitate accurate depth placement
- * Sliding stopper
- * Wire is marked to better control its placement
- * "V" shaped hook, 21G
- * Sterile, single-use
- * Latex-free



Item	Description	Pckg	Price
214.0001	Breast Localization Hook-Wire Needle - 21 Ga x 6 cm	25/box	\$419
214.0002	Breast Localization Hook-Wire Needle - 21 Ga x 9 cm	25/box	\$419
214.0003	Breast Localization Hook-Wire Needle - 21 Ga x 11 cm	25/box	\$419
214.0004	Breast Localization Hook-Wire Needle - 21 Ga x 15 cm	25/box	\$419

Radiation Protection Aprons

HOW TO CHOOSE YOUR APRON?

1. Choose the Style

Item	Description	Style
701.XXXX	Back Relief Apron	Front Protection
702.XXX	Buckle Apron	Front Protection
703.XXXX	Flexback Apron	Front Protection
704.XXXX	Surgical Drop Apron	Front Protection
705.XXXX	Tie Apron	Front Protection
706.XXXX	Velcro Adjustable Apron	Front Protection
707.XXXX	Vest Skirt Apron	Full Protection
708.XXXX	Vest Skirt Full Overlap Apron	Full Protection
709.XXXX	Wrap-Around Back Relief Apron	Full Protection
710.XXXX	Wrap-Around Back Relief FO Apron	Full Protection
711.XXXX	Wrap-Around Flexback Apron	Full Protection



2. Choose your Size

Gender	Size
Male	XS
Male	Small
Male	Medium
Male	Large
Male	XL
Male	2XL
Male	3XL
Male	4XL
Female	XS
Female	Small
Female	Medium
Female	Large
Female	XL
Female	2XL
Female	3XL
Female	4XL

3. Choose your Protection

Core Material	Lead Equivalence
Bilayer / Premier	0.25 Pb
Bilayer / Premier	0.35 Pb
Bilayer / Premier	0.50 Pb
Lead-Free	0.25 Pb
Lead-Free	0.35 Pb
Lead-Free	0.50 Pb
Ultra-Lite Lead	0.25 Pb
Ultra-Lite Lead	0.35 Pb
Ultra-Lite Lead	0.50 Pb
Lead	0.25 Pb
Lead	0.35 Pb
Lead	0.50 Pb



DOWNLOAD THE
SIZING CHART HERE

5. Add Personalization

- * Pocket
- * Embroidery & embroidery color
- * Sleeves
- * Logo
- * Binding Color

FIND MORE COLORS ONLINE!

4. Choose your Color


NYLON




RIPSTOP



VINYL



PRINT



Radiation Protection - Front Protection

Flexback Apron

- * Flex Elastic allows the apron's closure mechanism to "flex" away the apron's burden by shifting weight from shoulders to hips
- * Backstrap provides a secure, comfortable fit by not allowing the apron to slip on your shoulders
- * Closure: Velcro
- * Front Protection
- * Recommended for prolonged procedures

Item	Description	Pckg	Price
703.XXXX	Flexback Apron	Ea	From \$245



Back Relief Apron

- * Back relief belt transfers the weight of the apron from the shoulders to the hips
- * Increased comfort and back support
- * Closure: Buckle (belt)
- * Padded belt
- * Front Protection

Item	Description	Pckg	Price
701.XXXX	Back Relief Apron	Ea	From \$245



Buckle Apron

- * Belt is adjustable to accommodate different sizes
- * Criss-cross back straps provide a secure, comfortable fit
- * Easy-release buckle allowing for quick removal
- * Closure: Buckle (belt)
- * Front Protection

Item	Description	Pckg	Price
702.XXXX	Buckle Apron	Ea	From \$230



Surgical Drop Apron

- * Maintain your sterile field with the quick release Surgical Drop Apron
- * Easy-release velcro closures allow the surgical drop to fall even when underneath your scrubs
- * Padded shoulders and a criss-cross back strap provide maximum comfort
- * Closure: Velcro
- * Front Protection

Item	Description	Pckg	Price
704.XXXX	Surgical Drop Apron	Ea	From \$230



Tie Apron

- * Affordable option for maximum radiation protection
- * Criss crossing tie straps ensure a secure and comfortable fit
- * Closure: Tie
- * Front Protection

Item	Description	Pckg	Price
705.XXXX	Tie Apron	Ea	From \$220



Velcro Adjustable Apron

- * Crisscrossing back straps provide back support while securing the apron around your waist
- * Wide back straps prevent apron from slipping around your shoulder and back
- * Front Protection
- * Front Protection

Item	Description	Pckg	Price
706.XXXX	Velcro Adjustable Apron	Ea	From \$230



Thyroid Collars

- * Designed to fit higher along the neck
- * Complete coverage of the thyroid against scatter radiation
- * Closure: Velcro or Buckle
- * Optional Attachments: Loose, Clip, Key Ring

Item	Description	Pckg	Price
712.XXXX	Standard Thyroid Collar	Ea	\$50
712.XXXX	Original Thyroid Collar	Ea	\$50



Head Protection

- * Similar to a scrub cap
- * Provides lead protection against exposure to ionizing radiation
- * Adjustable velcro tab

Item	Description	Pckg	Price
713.XXXX	Head Cap	EA	From \$110
713.XXXX	ProGuard Smart Cap	Ea	From \$135
713.XXXX	Page Cap	Ea	From \$134



Radiation Protection - Full Protection

Vest & Skirt Apron

- * Two-piece design reduces back fatigue
- * Skirt comes with a premium non-slip fabric in the back
- * The vest features two overlapping panels for increased protection
- * Multiple buckle closures located on the vest & skirt pieces which provide additional adjustability for a more secure fit
- * Closure: Velcro & Buckles
- * Full Protection



Item	Description	Pckg	Price
707.XXX	Vest & Skirt Apron	Ea	From \$365

Vest & Skirt Full Overlap Apron

- * Vest features two 0.025mm Pb panels that create a 0.50mm overlap
- * Anti-Slip Fabric located in the back of the skirt
- * Multiple buckle closures located on the vest & skirt pieces which provide additional adjustability for a more secure fit
- * Closure: Velcro & Buckles
- * Full Protection



Item	Description	Pckg	Price
708.XXXX	Vest & Skirt Full Overlap Apron	Ea	From \$355

Wrap Around Back Relief Apron

- * Transfers the weight from the shoulders to the hips with the aid of a 6" bonded elastic belt
- * Double padded shoulders offer added comfort
- * Velcro break-away side seam for freedom of movement
- * Closure: Velcro (Standard) / Belt (Optional)
- * Full Protection



Item	Description	Pckg	Price
709.XXXX	Wrap-Around Back Relief Apron	Ea	From \$365

Wrap Around Full Overlap Apron

- * Transfers the weight from the shoulders to the hips with the aid of a 6" bonded elastic belt
- * Double padded shoulders offer added comfort
- * Wrap around eliminates redundant fabric and lead, reducing unnecessary weight
- * Closure: Velcro & Buckles
- * Full Protection



Item	Description	Pckg	Price
710.XXXX	Wrap-Around FO Apron	Ea	From \$365

Wrap Around Flexback Apron

- * Combines the best qualities of the flaxback apron with the full protection and back support of the wrap-around
- * Elastic stretch-back shifts weight from shoulders to hips
- * Tri-tab velcro closure for easy-on, easy-off solution
- * Full Protection



Item	Description	Pckg	Price
711.XXXX	Wrap-Around Flexback Apron	Ea	From \$405

Half Apron

- * Lead equivalence of 0.50 Pb
- * Small: 10" x 10"
- * Medium: 18" x 16"
- * Large: 24" x 20"
- * X-Large: 24" x 24"



Item	Description	Pckg	Price
714.XXXX	Half Apron	Ea	From \$70

Maternity Apron

- * Features easy ties and a quick release buckle to secure the back and front panels
- * Allows for complete adjustability during pregnancy
- * Includes an additional .50mm equivalent half apron for a total of 1.00mm Pb protection over the abdominal area
- * Full Protection
- * Closure: Velcro (Standard) / Buckle (Optional)



Item	Description	Pckg	Price
715.XXXX	Maternity Apron	Ea	From \$611

Gonad Shield with Belt

- * Comes with a one size fits all shield
- * Adjustable buckle also allows the user to fit almost any waist size while providing a comfortable fit
- * Lead equivalence of 0.50 Pb
- Closure - Buckle



Item	Description	Pckg	Price
716.XXXX	Gonad Shield	Ea	From \$45

Radiation Protection Glasses

9935 Ultralite

- * Lead equivalence 0.75 Pb
- * Weight: 60g
- * Available in Plano, Single Vision, Bifocal, and Progressive
- * Prescription Range: 5.00 to -5.00
- * Soft rubber nosebridge providing extra support and increased comfort
- * Available colors: Silver, Black, Red, Blue



Item	Description	Pckg	Price
750.XXXX	9935 Ultralite	Ea	From \$250

53 Wrap

- * Lead equivalence 0.75 Pb
- * Weight: 84g
- * Available in Plano, Single Vision, Bifocal, and Progressive
- * Prescription Range: 4.00 to -7.00
- * Side protection
- * Available colors: Red, Black, Tortoise



Item	Description	Pckg	Price
752.XXXX	53 Wrap	Ea	From \$160

70 Astroflex

- * Lead equivalence 0.75 Pb
- * Weight: 85g
- * Available in Plano, Single Vision, Bifocal, and Progressive
- * Prescription Range: 4.00 to -7.00
- * Side protection



Item	Description	Pckg	Price
755.XXXX	70 Astroflex w/ side shields	Ea	From \$145

Nike Brazen

- * Lead equivalence 0.75 Pb
- * Weight: 75g
- * Available in Plano, Single Vision, Bifocal, and Progressive
- * Prescription Range: 3.50 to -3.50
- * Available colors: Shiny Black, Matte Black, Matte Black/Blue, Tortoise



Item	Description	Pckg	Price
757.XXXX	Nike Brazen	Ea	From \$310

Ray Ban WayFarer

- * Lead equivalence 0.75 Pb
- * Weight: 62g
- * Available in Plano, Single Vision, Bifocal, and Progressive
- * Prescription Range: 3.00 to -3.00
- * Available colors: Matte Black



Item	Description	Pckg	Price
762.XXXX	Ray-Ban Wayfarer	Ea	From \$346

9941 Ultralite

- * Lead equivalence 0.75 Pb
- * Weight: 64g
- * Available in Plano, Single Vision, Bifocal, and Progressive
- * Prescription Range: 6.00 to -9.00
- * Curved front and temples with rubber detail assures a close secure fit
- * Available colors: Silver, Black, Red, Blue



Item	Description	Pckg	Price
751.XXXX	9941 Ultralite	Ea	From \$250

42 Fitover

- * Lead equivalence 0.75 Pb
- * Made to wear over prescription eyewear
- * Weight: 90g
- * Available in Plano only
- * Side protection



Item	Description	Pckg	Price
753.0001	42 Fitover w/ Lead Strip Side Shield	Ea	\$180

ProGuard Leader

- * Soft ventilated nose bridge with rubberized temples
- * Lead equivalence 0.75 Pb
- * Weight: 70g
- * Available in Plano, Single Vision, Bifocal, and Progressive
- * Prescription Range: 3.50 to -3.50
- * Side protection 0.50 Pb



Item	Description	Pckg	Price
756.XXXX	ProGuard Leaders	Ea	From \$250

Oakley Crosslink

- * Lead equivalence 0.75 Pb
- * Weight: 62g
- * Available in Plano, Single Vision, Bifocal, and Progressive
- * Prescription Range: 3.50 to -9.00
- * Available colors: Satin Black, Matte Grey



Item	Description	Pckg	Price
758.XXXX	Oakley Crosslink	Ea	From \$376


Ray Ban LightRay

- * Lead equivalence 0.75 Pb
- * Weight: 46g
- * Available in Plano, Single Vision, Bifocal, and Progressive
- * Prescription Range: 3.00 to -5.00
- * Available colors: Matte Black, Aqua Blue



Item	Description	Pckg	Price
763.XXXX	Ray-Ban LightRay	Ea	From \$375

Distributed By:



*Not available for ultrasound gels, lubricants, and disinfectants

Phone (877) 220-4567 - Fax (786) 231-0336
www.edm-us.com - info@edm-us.com